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
Figure 1. Real-time PCR reproducibility. This plot shows results from amplification of KAZ plasmid DNA in 10-fold dilutions (A to J: starting with A = 10^9 to J = 10^0 , no template control (NTC)) using the 384-well block and Applied Biosystems™ TaqMan® Fast Advanced Master Mix. The data show highly reproducible results over 10 logarithmic units of input template amount, illustrating the broad linear dynamic range of the system.

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
Figure 2. Comparison of $\Delta\Delta C_t$ data between QuantStudio 7 Pro, QuantStudio 7 Flex, and 7900HT systems demonstrates equivalent relative expression level of GFAP, ALB, and KRT1 in cDNA of universal human reference RNA, total liver RNA (LIV), and total brain (BRN) RNA.


Figure 3. Primary analysis of a completed run. The amplification curves, assay maps, and results including quality control (QC) alerts are displayed after a run is complete. The curves can be adjusted as needed, similar to other QuantStudio instruments.

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
Trade in larger and older instruments ...

... and save lab bench space


54.7 x 33.8 x 52.5 cm
(H x W x D), and the screen is in full color.

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

Find out more at thermofisher.com/taqman-rfid

Description	No. of assays	No. of samples	Cat. No. (96-well 0.1 mL)	Cat. No. (96-well 0.2 mL)
TaqMan Array Plates with RFID	N/A	N/A	A43823	
TaqMan Array, Standard, 96-well plate; Format 8	8	12	4413263	4413266
TaqMan Array, Standard, 96-well plate; Format 16/16 Plus	16	6	4413261/4413262	4413264/4413265
TaqMan Array, Standard, 96-well plate; Format 32/32 Plus	32	3	4413259/4413260	4391528/4391529
TaqMan Array, Standard, 96-well plate; Format 48/48 Plus	48	2	4413257/4413258	4391526/4391527
TaqMan Array, Standard, 96-well plate; Format 96/96 Plus	96	1	4413255/4413256	4391524/4391525

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- Remote monitoring that shows run status and amplification plots in real time
- Data sharing and analysis features, using the cloud-enabled Connect platform
- Planning and scheduling your next experiment without delay for maximum productivity
- Be the ultimate qPCR pro in your lab or anywhere

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	QuantStudio 6 Pro Real-Time PCR System	QuantStudio 7 Pro Real-Time PCR System
样品容量 (wells)	96 or 384	
反应体积*	96-well 0.1 mL block: 10–30 μ L 96-well 0.2 mL block: 10–100 μ L 384-well block: 5–20 μ L	96-well 0.1 mL block: 10–30 μ L 96-well 0.2 mL block: 10–100 μ L 384-well block: 5–20 μ L TaqMan Array Card: ~1.5 μ L
激发光源	Bright white LED	
滤光片或颜色组合数	5	21
激发光波长/发射光波长	450–680 nm/500–730 nm	
多重荧光定量靶标种数 (Multiplex qPCR)	5 or 6 targets	
反应模块最快变温速率	96-well 0.1 mL block: 9.0°C/sec 96-well 0.2 mL block: 6.5°C/sec 384-well block: 6.0°C/sec	96-well 0.1 mL block: 9.0°C/sec 96-well 0.2 mL block: 6.5°C/sec 384-well block: 6.0°C/sec TaqMan Array Card: 3.0°C/sec
样品平均变温速率	3.66°C/sec	
模块温度均一性	0.4°C	
温控范围	4–99.9°C	
所支持的燃料类型	FAM™, SYBR™ Green, VIC™, JOE™, TET™, ABY™, NED™, TAMRA™, Cy®3, JUN™, ROX™, Texas Red™, Mustang Purple™, Cy®5, LIZ™, Cy®5.5	
数据电子签名功能 (SAE)	No	Yes, with no additional fees
自动机械臂整合功能	No	Yes

* Blocks are being released at intervals; some may not be currently available.

For more information and full testimonials, go to thermofisher.com/quantstudiopro

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
Figure 4. Comparative data as shown with onboard software. This chart shows duplex comparative C_t data with results from the Design and Analysis Software. Experiment done with a reaction volume of 10 μ L using RNase P, TaqMan Fast Advanced Master Mix, and Applied Biosystems™ Endura™ 384-well Optical Plate.


Figure 5. Relative quantification. These charts enable efficient data review in the Design and Analysis Software.

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人性化—automatically load your settings and plate setup; no passwords required to log in, thanks to facial authentication

高效率—facilitate hands-free operations with voice commands; maximize uptime with Smart Help feature to access technical support experts for more efficient troubleshooting and resolution; reduce downtime using Smart Remote Support collaboration tool to guide you through repairs remotely; eliminate manual steps to obtain plate layout, protocol, and assay information on the instrument with TaqMan Array Plates with RFID

简单易用—straightforward block changes with no tools required; streamlined workflow directly from the touchscreen

富有成效—built-in help videos for commonly used activities; increased ergonomics with larger touchscreen, improved screen angle and optional Bluetooth™ keyboard; access data anytime and anywhere with cloud-enabled services from Connect

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Find out more at thermofisher.com/instrumentservices

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* Extended warranty packages include the instrument, SmartStart Orientation, and a 1-year AB Assurance Plan with 1 planned maintenance (PM) visit. Packages are not available in all countries. Contact your local sales representative for more information.

** Some block configurations may not yet be available. Go to thermofisher.com/quantstudiopro for a complete list and availability.

Ordering information

Product	Size	Cat. No.
Consumables		
TaqMan Fast Advanced Master Mix	5 mL	4444557
	50 mL	4444558
PowerTrack SYBR Green Master Mix	5 mL	A46109
	50 mL	A46113
High Capacity RNA-to-cDNA Kit	50 reactions	4387406
Invitrogen SuperScript IV VIL0 Master Mix	500 reactions	11756500
Instrument Qualification (IQ), Operational Qualification (OQ), and Instrument Performance Verification (IPV) services*		
IQ/OQ/IPV service for QuantStudio 6 or 7 Pro systems, 96-well, 0.1 mL block		A45610
IQ/OQ/IPV service for QuantStudio 6 or 7 Pro systems, 96-well, 0.2 mL block		A44230
IQ/OQ/IPV service for QuantStudio 6 or 7 Pro systems, 384-well block		A45576
IQ/OQ/IPV service for QuantStudio 7 Pro systems, TAC block		A47663

* OQ/IPV service is not available in all countries. Contact your local sales representative for availability information.

Experience faster, easier qPCR in a compact design. See the future of qPCR at thermofisher.com/quantstudiopro

ThermoFisher
SCIENTIFIC