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SpectraMax iD5 Multi-Mode Microplate Reader

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
*View data on the large touchscreen for quick reading of results*

*Export data to a USB drive for analysis in the program of your choice*

*Automatically send data to any workstation on the same network*

*Use SoftMax Pro 7 Software for advanced data acquisition and analysis*

# Reliable performance with a proven track record

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


An ATP standard curve spanning five decades was run using the ATPlite 1step Luminescence Assay System (PerkinElmer) on the SpectraMax iD5 reader. Standard concentrations ranged from  $1 \times 10^{-11}$  M to  $1 \times 10^{-6}$  M. A wide linear dynamic range ensures accurate assay results across a broad span of sample types. Standards were plotted using a log-log curve fit in SoftMax Pro Software.

The SpectraMax iD5 reader is fully compatible with the SpectraDrop™ Micro-Volume Microplate, enabling quantitation of precious low-volume samples. 4- $\mu$ L DNA standards from 2 ng/ $\mu$ L to 1000 ng/ $\mu$ L were read in absorbance detection mode with a preconfigured protocol in SoftMax Pro Software. Performance matching the 2 ng/ $\mu$ L sensitivity specification is demonstrated here.

## 技术规格 – SpectraMax iD5 reader

### General specifications

Dimensions (in.)	15.79 (H) x 20.94 (W) x 23.54 (D)
Dimensions (cm)	40.1 (H) x 53.2 (W) x 59.8 (D)
Weight	88.1 lbs. (40 kg)
Power requirements	100–240 VAC, 2 A, 50/60 Hz
Robotic compatible	Yes

### General performance

Plate formats	6 to 384 wells
Light source	Xenon flash lamp
Reading capabilities	Microplates, cuvettes (via adapter)
Detectors	Photomultiplier Tube and Photodiode
Shaking	Linear, orbital and double orbital
Temp. control	5°C above ambient to 66°C <sup>①</sup>
Temp. uniformity	± 0.75°C
Temp. accuracy	± 1°C at 37°C set point
Spectral scanning	Abs, FI, Lum, TRF
Endpoint reading	Abs, FI, Lum, TRF, FP, FRET, TR-FRET
Kinetic reading	Abs, FI, Lum, TRF, FP, FRET, TR-FRET
Well scanning	Over 20 by 20
Wavelength selection	1.0 nm increments

### Standard read times (minutes:seconds)

	96 wells	384 wells
Absorbance	0:30	1:30
Fluorescence intensity <sup>②</sup>	0:30	1:30
Luminescence <sup>③</sup>	0:30	1:30
Time-Resolved Fluorescence	0:30	1:30
Fluorescence Polarization	1:00	3:00

### Absorbance photometric performance

Wavelength range	230–1000 nm
Wavelength bandwidth	4.0 nm full width half maximum
Wavelength accuracy	± 2.0 nm across wavelength range
Wavelength repeatability	± 1.0 nm
Photometric range	0–4.0 OD
Photometric resolution	0.001 OD
Photometric accuracy	< ±0.010 OD ±1.0% 0–3 OD VIS, 0–3 OD UV
Photometric precision (repeatability)	< ±0.003 OD ±1.0% 0–3 OD VIS, 0–3 OD UV
Stray light	< 0.05% @ 260 nm, 280 nm

### Fluorescence intensity performance

Wavelength range (EX mono)	250–830 nm
Wavelength range (EM mono)	270–850 nm
Wavelength range (Filters)	See accessory list
Wavelength selection (mono)	1.0 nm increments
Dynamic range	> 6 logs

### Optimized top sensitivity (fluorescein)

	Monochromator	Filters
96 wells	1 pM	0.3 pM
384 wells	1 pM	0.5 pM

### Optimized bottom sensitivity (fluorescein)

	Monochromator	Filters
96 wells	2 pM	2 pM
384 wells	2.5 pM	2.5 pM

### Fluorescence polarization performance

Wavelength range (EX mono)	(300-)④ 400–750 nm, 1.0 nm increments
Wavelength range (EM mono)	(300-)④ 400–750 nm, 1.0 nm increments
Wavelength range (Filters)	See accessory list

### Detection limit<sup>⑤</sup>

	Monochromator	Filters
96 wells black	2 mP @ 10 nM	1 mP @ 1 nM
384 wells black	2 mP @ 10 nM	2 mP @ 1 nM
Measurement range <sup>⑥</sup>	Delta > 200 mP	Delta > 320 mP

### Luminescence performance

Wavelength range	300–850 nm 300–650 nm for “All Wavelengths” setting
Wavelength selection	Choice of simultaneous detection of all wavelengths or selection in 1.0 nm increments
Dynamic range	> 7 decades
Cross-talk	< 0.1% in white 96- and < 0.2% in 384-well microplates
Detection limit	20 amol ATP (“Flash” luminescence using Promega ENLITEN® ATP Assay System)

### Time-Resolved Fluorescence performance

Wavelength range (EM mono)	450–750 nm
Wavelength range (EX filter)	350 nm
Wavelength range (EM filter)	490 nm (Terbium), 616 nm (Europium) For other filters, please see accessory list
Linear dynamic range	Up to 5 logs

### Detection limit

	96 wells (white)	384 wells (white)
Standard TRF	30 fM Europium (6 amol/well)	30 fM Europium (3 amol/well)
Enhanced TRF Module (Optional)	10 fM Europium (2 amol/well)	10 fM Europium (1 amol/well)

### Injector system with SmartInject Technology (optional)

Injectors	2
Read modes	Absorbance, fluorescence, luminescence
Dispense accuracy	± 5% at 100 µL
Dispense precision	CV ≤ 2% at 100 µL
Dead volume	Injector Tubing: 250 µL < 10 µL with Reverse Prime function

<sup>①</sup> For > 66°C, minimum 25°C ambient temperature is required.

<sup>②</sup> 10 msec integration time for fluorescence measurement. 100 msec integration time for luminescence for a 96-well plate and 40 msec integration time for a 384-well plate.

<sup>③</sup> Requires optional UVIS polarizer

<sup>④</sup> 1x Stdev Fluorescein replicates [mP]

<sup>⑤</sup> ThermoFisher P3088, FP One-Step reference kit