QuantStudio 7 Flex 和 QuantStudio 6 Flex

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See our growing catalog of sensitive, specific real-time PCR probe and primer sets for human viruses and other areas of infectious disease research.

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Detect multiple bacteria in the same run, including *Salmonella*, *Campylobacter*, *E. coli* O157:H7, *Listeria monocytogenes*.

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Designed to detect and monitor waterborne pathogens in recreational and drinking water supplies.

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Designed to detect mycoplasmas, viruses, and residual host cell contamination for pharmaceutical, cosmetics, and personal care product manufacturing.

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A growing menu of new and valuable content for you to use on your high-throughput real-time PCR systems, including pharmacogenomics, vaginal microbiota, and *CFTR* mutation analysis.

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Solutions for analyzing stem cells, determining stemness, and studying gene regulation and translation in stem cells.

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Pre-designed Applied Biosystems™ TaqMan® Assays for more than 175 ADME and CYP targets, including >95% of ADME core markers and a warfarin metabolism panel.

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Enabling robust, reliable detection and quantitation of markers for cancer and genetic diseases.

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Instruments, reagents, and kits designed for plant researchers that enable remarkable agricultural discoveries—from improved crops that feed more people to sustainable biofuels.

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	QuantStudio 6 Flex system	QuantStudio 7 Flex system
Block configuration	96-well, 96-well Fast, 384-well	96-well, 96-well Fast, 384-well, TaqMan Array Card
Dimensions (H x W x D)	73.8 x 50.1 x 66.1 cm	73.8 x 50.1 x 66.1 cm
Block change design	Block change from front in less than 1 min; no tools required	Block change from front in less than 1 min; no tools required
Run time	30 min (96-well Fast) 35 min (384-well)	30 min (96-well Fast) 35 min (384-well)
Ramp-up rates	3.9°C/sec (96-well) 6.5°C/sec (96-well Fast) 3.5°C/sec (384-well)	3.9°C/sec (96-well) 6.5°C/sec (96-well Fast) 3.5°C/sec (384-well) 4.3°C/sec (TaqMan Array Card)
Ramp-down rates	3.6°C/sec (96-well) 6.0°C/sec (96-well Fast) 3.0°C/sec (384-well)	3.6°C/sec (96-well) 6.0°C/sec (96-well Fast) 3.0°C/sec (384-well) 3.7°C/sec (TaqMan Array Card)
Well-to-well variability	±0.25°C	±0.25°C
Excitation source	OptiFlex System lamp	OptiFlex System lamp
Detection channels	Excitation 455–672 nm Emission 505–723 nm	Excitation 455–672 nm Emission 505–723 nm

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Product	Cat. No.	Instrument + 1-year extended warranty with AB Assurance Cat. No.*
QuantStudio 6 Flex 96-Well Real-Time PCR System, laptop configuration	4485689	A27163
QuantStudio 6 Flex 96-Well Fast Real-Time PCR System, laptop configuration	4485699	A27161
QuantStudio 6 Flex 384-Well Real-Time PCR System, laptop configuration	4485691	A27960
QuantStudio 6 Flex 96-Well Real-Time PCR System, desktop configuration	4485692	A27146
QuantStudio 6 Flex 96-Well Fast Real-Time PCR System, desktop configuration	4485697	A27162
QuantStudio 6 Flex 384-Well Real-Time PCR System, desktop configuration	4485694	A27958
QuantStudio 7 Flex 96-Well Real-Time PCR System, laptop configuration	4485688	A27154
QuantStudio 7 Flex 96-Well Fast Real-Time PCR System, laptop configuration	4485698	A27155
QuantStudio 7 Flex 384-Well Real-Time PCR System, laptop configuration	4485695	A27152
QuantStudio 7 Flex TaqMan Array Card Real-Time PCR System, laptop configuration	4485700	A27150
QuantStudio 7 Flex 96-Well Real-Time PCR System, desktop configuration	4485690	A27157
QuantStudio 7 Flex 96-Well Fast Real-Time PCR System, desktop configuration	4485693	A27156
QuantStudio 7 Flex 384-Well Real-Time PCR System, desktop configuration	4485701	A27153
QuantStudio 7 Flex TaqMan Array Card Real-Time PCR System, desktop configuration	4485696	A27151
QuantStudio 6/7 Flex 96-Well Block Upgrade Kit	4453543	NA
QuantStudio 6/7 Flex 96-Well Fast Block Upgrade Kit	4453544	NA
QuantStudio 6/7 Flex 384-Well Block Upgrade Kit	4453545	NA
QuantStudio 7 Flex TaqMan Array Card Block Upgrade Kit	4453546	NA

* Includes SmartStart orientation.

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96-well block	Cat. No.	96-well block Fast	Cat. No.	384-well block	Cat. No.
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plates with Barcode	4483354	MicroAmp EnduraPlate Optical 96-Well Fast Clear Reaction Plates with Barcode	4483485	MicroAmp EnduraPlate Optical 384-Well Clear Reaction Plates with Barcode	4483285
MicroAmp Optical 96-Well Reaction Plate with Barcode	4306737	MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL	4346906	MicroAmp Optical 384-Well Reaction Plate with Barcode	4309849
MicroAmp Optical 8-Tube Strip with Attached Optical Caps, 0.2 mL	A30588	MicroAmp Fast 8-Tube Strip	4358293	MicroAmp Optical 384-Well Reaction Plate	4343370
MicroAmp Optical Adhesive Film	4360954	MicroAmp Optical 8-Cap Strips	4323032	MicroAmp Optical Adhesive Film	4360954
		MicroAmp Optical Adhesive Film	4360954		

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	维修中心服务计划*			现场服务计划	
	AB Repair Center Support Plus	AB Repair Center Support Plus Care	AB Maintenance Plus	AB Assurance	AB Complete
On-site response time			Target 2 business days**	Guaranteed 2 business days**	Guaranteed next business day**
Scheduled on-site planned maintenance (PM)		✓	✓	✓	✓
Remote diagnostics	✓	✓	✓	✓	✓
Parts, labor, and travel for repair included	✓	✓	10% discount optional add-on in selected regions	✓	✓
Computer repair and replacement included	✓	✓		✓	✓
Priority access to Tech Support (Mon–Fri, 8 a.m.—5 p.m. local time)	✓	✓		✓	✓
Priority access to Remote Service Engineer	✓	✓		✓	✓
Requalification post-PM and critical repairs					✓
Field Application Scientist consultation					✓
Loaner instrument issued during repair (Repair Center plans only)	✓	✓			

* Repair Center service plans available for QuantStudio 3 and 5 Real-Time PCR Systems only.

** Response times vary by region.